

# Radnoti

## Radnoti 24mm Heart Chamber and Buffer Reservoir Adapters:



140131



140132



The Radnoti #120147 adapter fits between the heart chamber and the buffer reservoir. The two Luer side ports permit the addition of new medium or the sampling of effluent, as well as permitting the gassing of the heart chamber with nitrogen during hypoxia experiments. Constructed of borosilicate glass. The Radnoti #140132 adapter converts the heart chambers 24mm ground glass exit port to a female Luer fitting for easy connection of tubing for the collection of or draining of chamber fluids. The Radnoti #140131 adapts the buffer reservoir Radnoti #120146 to a female luer fitting for easy connection of tubing for the collection

Catalog #	Description	Pkg.
120147	24mm heart chamber adapter	1
140131	24mm buffer reservoir adapter	1
140132	24mm male to female adapter	1



120147

### Radnoti Post Heart Buffer Reservoir

The Radnoti #120147 adapter fits between the heart chamber and the buffer reservoir. The two Luer side ports permit the addition of new medium or the sampling of effluent, as well as permitting the gassing of the heart chamber with nitrogen during hypoxia experiments. Constructed of borosilicate glass. The Radnoti #140132 adapter converts the heart chambers 24mm ground glass exit port to a female Luer fitting for easy connection of tubing for the collection of or draining of chamber fluids. The Radnoti #140131 adapts the buffer reservoir Radnoti #120146 to a female luer fitting for easy connection of tubing for the collection

Catalog #	Description	Pkg.
120146	Post Heart Buffer Reservoir	1



### Radnoti Glass Technology, Inc

227 W. Maple Ave. Monrovia, CA 91016  
 (800) 428-1416 (626) 357-8827 fax (626) 303-2998  
[www.radnoti.com](http://www.radnoti.com)